

SpecClean™

**ULTRASONIC CLEANER
SOLUTION**

Neutral pH



Concentrated Ultrasonic Cleaner Solution

Spec Clean™ Ultrasonic Cleaner is a super concentrated instrument cleaning solution formulated for use in ultrasonic cleaning devices. The biodegradable solution meets the demand for more environmentally friendly products.

Intended Use

Spec Clean Ultrasonic Cleaner is intended for cleaning instruments by aiding the cavitation process of ultrasonic cleaners. It is designed to minimize manual handling/scrubbing of contaminated items and is especially useful for cleaning instruments with serrated edges, joints and hinges. It will safely clean delicate instruments, anodized aluminum and carbon steel curets or burs.

Enhances Ultrasonic Cleaning Action

- Non-foaming formula best for ultrasonic cleaners.
- Penetrating agents improve cleaning action in small crevices and hinges.
- Odor suppressant neutralizes organic and chemical odors.
- Economical super concentrated 1 oz. unit dose packet may be diluted with ordinary tap water to make one gallon of cleaning solution.



Helping improve the health and quality of life for animals.

Protects Delicate Instruments:

- Anti-corrosive agents protects metal surfaces.
- Best results for carbon steel items including burs and curets.
- Safe cleaning for jointed orthodontic pliers.
- Enhances instrument appearance and shine.

Meet CDC Guidelines for Instrument Decontamination

- Effective cleaning is essential prior to sterilization.
- Effective ultrasonic cleaning reduces risk of injury during manual scrubbing of instruments.

Product Capabilities

Spec Clean Ultrasonic Cleaner Solution rapidly removes organic and inorganic soils. The detergent formula lifts and disperses soils during the ultrasonic cleaning process. A powerful anti-corrosive system protects delicate instruments including carbon steel and anodized aluminum. The anti-corrosive agents are especially effective to protect box locks in jointed items and hinges.

Use Precautions

Do not add other chemicals to Spec Clean Ultrasonic Cleaner Solution. Change the cleaning solution daily or more frequently if soiled or cloudy. Empty the ultrasonic tank, rinse and wipe out completely DAILY. Do not allow residue to accumulate. Follow ultrasonic equipment manufacturer's recommendations for cleaning and maintenance.

Due to the corrosive effects of water, it is not recommended to soak metal items for prolonged periods (over two hours or overnight).

Spec Clean Ultrasonic Cleaner is an eye irritant and mild skin irritant. Follow CDC Guidelines for protective attire during instrument decontamination: safety goggles and/or face shield, mask (for biological hazards) and water resistant gloves. Store in cool, dry area away from direct sunlight. Solution will fade when exposed to direct sunlight.

Availability

Stocked by MicroCare Medical distributors. Contact MicroCare Medical for a distributor near you or visit MicroCare.com

Product Specifications

Dilution Ratio	1 oz. or 1 unit dose per gallon
pH Concentrate	7.0-8.0
Temp. of use	68°F (20°C) to 140°F (60°C). Do not exceed 140°F (60°C). Avoid running continuous cycles that may overheat the cleaning solution
Odor	Mint
Shelf life	3 years from date of manufacture
Disposal	Dispose according to all local, state and federal regulations. This product contains no phosphates

Package Sizes

Size	Case Qty	Order #
1 Gallon	4	VGVUSNP128-1

Rev. VGVUS/TDS 23319

© 2023 MicroCare, LLC. All Rights Reserved. "MicroCare Medical", "Spec Clean" are trademarks or registered trademarks of MicroCare, LLC. The information set forth herein is based on data believed to be reliable. MicroCare Medical makes no warranties express or implied as to its accuracy and assumes no liability arising out of its use by others. This publication is not to be taken as a license to operate under, nor infringe upon, any patents not herein expressly described.



6120 East 58th Avenue, Commerce City, CO 80022
(800) 843-3343 • www.microcare.com